

# Feasibility Evaluation of Blow Fill Seal Process and Compatibility with Aluminum Phosphate Adjuvanted Recombinant RSV F Nanoparticle Vaccine

YH Lin<sup>1</sup>, E Orvisky<sup>1</sup>, R Hau<sup>1</sup>, S Boddapati<sup>1</sup>, K Patel<sup>1</sup>,  
E Wu<sup>1</sup>, S Gschwind<sup>3</sup>, Y Schwander<sup>3</sup>, T Kram<sup>2</sup>, T Hahn<sup>1</sup>  
<sup>1</sup>Novavax, Inc., Gaithersburg, MD; <sup>2</sup>Rommelag USA Inc., <sup>3</sup>Maropack AG

Contact:  
Yen-Huei Lin, Ph.D.  
ylin@Novavax.com  
Phone: 240.268.2000



## INTRODUCTION

Respiratory syncytial virus (RSV) is one of the common causes of childhood acute lower respiratory tract infections (ALRI) in infants and young children worldwide.

In a 2015 analysis, six million estimated episodes of RSV-severe ALRI occurred in children younger than 5 years of age in low-income and middle-income countries (LMIC). The overall RSV-ALRI mortality was estimated to be as high as 118,200. In children younger than six months, 1.4 million hospital admissions and 27,300 in-hospital deaths were due to RSV-ALRI<sup>1</sup>.

ResVax (Pre-fusogenic RSV F nanoparticle vaccine) is currently being assessed in the Prepare™ (Phase 3) trial for the protection of infants via maternal immunization in healthy third trimester pregnant women.

To ensure cost effective and cold chain efficient vaccine for LMIC, a novel product presentation is needed in contrast to conventional glass vials and syringes.




## CONCLUSIONS

- Product uniformity can be achieved using the BFS technology.
- The stability profile of RSV F nanoparticle vaccine in BFS under accelerated and intended storage conditions appears comparable to the profile in glass vials and syringes.
- The result of this feasibility assessment provides data to enable a potential cost-effective product presentation of the RSV F nanoparticle vaccine for WHO pre-qualification pathway.

## REFERENCES

- Shi T *et al.* Global, regional, and national disease burden estimates of acute lower respiratory infections due to respiratory syncytial virus in young children in 2015: a systematic review and modelling study. *Lancet* 2017 390(10098): 946-958.
- PATH The global respiratory syncytial virus burden 2017. <https://www.path.org/publications/detail.php?i=2810>; accessed Feb 2018
- PDA Technology Report No. 77: The Manufacture of Sterile Pharmaceutical Products Using Blow-Fill-Seal Technology 2017.


## BLOW-FILL-SEAL (BFS) TECHNOLOGY




| Estimates  | 2mL Vial                      | 1mL Standard Pre-filled syringe | Blow-Fill-Seal containers     |
|--|-------------------------------|---------------------------------|-------------------------------|
| Relative Fill Process & Components Cost          | Medium                        | High                            | Low                           |
| Relative Product loss (process, container & use) | High                          | Low                             | Low                           |
| Primary container weight                         | 6-7 g                         | 6-7 g                           | 1 g                           |
| Cold Chain Volume                                | ~20 cm <sup>3</sup> /10 doses | ~60 cm <sup>3</sup> /10 doses   | ~10 cm <sup>3</sup> /10 doses |

- BFS technology<sup>3</sup> has been a robust method to produce aseptic pharmaceuticals since the late 1960's.
- In the past, the technology was not generally feasible for protein based biopharmaceuticals and vaccines due to elevated process temperature.
- Drivers to use BFS technology for vaccines include high throughput production, single inventory supply chain (only resin; no vial or syringe components), and a lower storage volume requirement.

### Advances in Controlling BFS Process Temperature Enable Application for Protein Based Biopharmaceuticals and Vaccines



|                             |                   |
|-----------------------------|-------------------|
| Formulation Temperature     | 8°C               |
| Filling Volume              | 0.34 ml...1.67 ml |
| Wall Thickness              | 0.4 mm...0.7 mm   |
| Mould & Formulation Cooling | on                |



The data above show the temperature profiles within the mould and inside the ampoule. The BFS process with cooling shows that temperature can be controlled within a range that is acceptable for proteins.

## ACKNOWLEDGEMENT

- Bill & Melinda Gates Foundation for grants to Rommelag and Maropack that supports this study.
- Ray Prasad, Tina Lorensen, and Torey de Rozario at Bill & Melinda Gates Foundation; Russel Wilson and Matt Lawlor at Novavax for their support and discussion to initiate this feasibility study.
- Matthew Peters and Nels Peterson at Global Good for their work on compact BFS container design and BFS fill recirculation system.
- Greg Glenn, Amy Fix, Lou Fries, Erika Trahan, and Robert Fuentes for their review and comments; Caitlin Cavazos for her assistance in preparing this poster.

## FEASIBILITY ASSESSMENT OF ALUMINUM PHOSPHATE ADJUVANTED RSV F VACCINE WITH BLOW-FILL-SEAL

### Objective


- Evaluate BFS technology for product compatibility and process feasibility with aluminum phosphate adjuvanted RSV F nanoparticle vaccine.
- Enable data driven decision for potential cost effective product presentation for WHO pre-qualification pathway.

### Process Control Challenge with Aluminum Phosphate Adjuvanted Product




To avoid product uniformity issues aluminum phosphate particles need to be continuously mixed in a product holding tank and continuously recirculated in the fill line.

### Process Control for Fill Uniformity with BFS Technology Is Feasible



### Fill Uniformity Further Improved with Recirculation Loop



### Stability Testing with and without Moisture/Gas Barrier Secondary Packaging


| Temp  | Test       | Time 0 | 1M | 3M | 6M | 9M | 12M | 18M | 24M |
|-------|------------|--------|----|----|----|----|-----|-----|-----|
| 2-8°C | ELISA      | x      | x  | x  | x  | x  | x   | x   | x   |
|       | A280       | x      | x  | x  | x  | x  | x   | x   | x   |
|       | Appearance | x      | x  | x  | x  | x  | x   | x   | x   |
|       | RP-HPLC    | x      | x  | x  | x  | x  | x   | x   | x   |
| 25°C  | ELISA      | x      | x  |    |    |    |     |     |     |
|       | A280       | x      | x  |    |    |    |     |     |     |
|       | Appearance | x      | x  |    |    |    |     |     |     |
|       | RP-HPLC    | x      | x  |    |    |    |     |     |     |
| 37°C  | ELISA      | x      |    |    |    |    |     |     |     |
|       | A280       | x      |    |    |    |    |     |     |     |
|       | Appearance | x      |    |    |    |    |     |     |     |
|       | RP-HPLC    | x      |    |    |    |    |     |     |     |



Limited water loss even if BFS containers were not stored in gas / moisture barrier secondary packaging

| Storage Conditions | Water Loss per Month |
|--------------------|----------------------|
| 2-8 °C             | 0.03%                |
| 25 °C              | 0.39%                |
| 37 °C              | 1.35%                |

### RSV F Nanoparticle Vaccine Stability at 2-8 °C



### RSV F Nanoparticle Vaccine Stability at 25 and 37 °C

